




GUNT experimental flumes Laboratory design

experimental flumes including the water tank.

GUNT will gladly undertake the precise laboratory planning for you to set up the experimental flumes.

The following table lists the space requirements of all GUNT A lifting device is recommended when placing larger models in the experimental sections of HM 161.



	А	B (excl. G)	С	C (incl. G)	D	Е	F	G	н	Height B (excl. H)	Height B (incl. H)	Required room height
HM 160	2,5 m 5,0 m	4,3 m 6,9 m	0,7m		1,0 m	1,5 m (>1 m)	2,0 m			1,35m	1,80m	2,3m
HM 162/ HM 163	5,0 m 7,5 m 10,0 m 12,5 m	9,2 m 11,7 m 13,6 m 16,0 m	1,0 m 1,0 m 2,2 m 2,2 m	2,2m 2,2m 2,2m 2,2m	1,0 m	1,5m (>1m)	2,5 m	1,0m	1,7m	2,20m	2,90m	with sediment feeder: min. 4,5 m
HM 161	16,0 m	22,0 m	4,0 m	4,0 m	2,0 m	1,5 m (>1 m)	1,0 m	1,0 m	in C incl.	2,70m	3,70m	with sediment feeder: min. 5 m

Installation requirements

This section provides some tips for planning a laboratory in which an experimental flume is going to be set up:

- the laboratory should be on the ground floor
- the floor must have sufficient load capacity
- the floor and the skirting area of the walls should be water-resistant
- lacktriangleright the transportation routes to and within the laboratory must be spacious
- the water supply and drains must be big enough for larger amounts of water
- the two larger experimental flumes HM 162, HM 163, and HM 161 require three-phase alternating current

An example of laboratory planning

The drawing below shows the planning for a laboratory that contains the experimental flume HM 162 (10 m long experimental section), a few other GUNT units on fluid mechanics and workstations for the students.

In this case the models for HM162 are stored on tables.

A small cabinet in the corner contains tools and can be used to store instruction manuals.

